







The CURIS® flow irrigation pump seamlessly integrates into surgical workflows. Moreover, the CURIS® flow is compatible with various irrigation tubing sets available on the market.



## **INCLUDED:**



1x 360901 CURIS® flow IV pole mounting bracket

1x 93008120 CURIS® RF Generator interconnecting cable

## **OPTIONAL ACCESSORIES:**



**360115** CURIS<sup>®</sup> flow one-pedal footswitch with push-button

Technical Data	
Mains connection	100 - 240 V; 50 / 60 Hz
Power consumption	Max. 30 VA
Protection class	1
Туре	CF (cardiac floating); defibrillator safe
Protection level	IP21
Weight	Approx. 3.0 kg
Dimensions	W x H x D 230 mm x 125 mm x 250 mm
Standards	IEC 60601-1 IEC 60601-1-2

Product availability is subject to regulatory approval in individual markets. Products may therefore not be available in all markets. The listed lengths and sizes serve as a guideline and may be rounded up or down. The actual lengths may vary slightly.



PRECISION ELECTROSURGERY Made in Germany SUTTER MEDICAL TECHNOLOGIES USA, INC. 3159 ROYAL DRIVE · SUITE 340 - ALPHARETTA, GA 30022 WWW.SUTTER-USA.COM · INFO@SUTTER-USA.COM