



Irrigating Bipolar Forceps



Irrigating Bipolar Forceps

Irrigating Bipolar Forceps · YASARGIL non-stick · irrigation · stop-pin

non-stick technology
Silver alloy tips prevent overheating of the tips during coagulation. The non-stick effect is permanently ensured, as the tips are made of solid noble metals.

Stop-pin
Minimizes scissoring of the tips.



Stop-pin

Irrigation cap
A universal hub connector allows irrigation tubing to be connected easily.

Handle
The ergonomically designed handle provides less fatigue and a well balanced instrument.



0.4 mm
SUT-450SP R

0.7 mm
SUT-453SP R

1.0 mm
SUT-455SP R

1.2 mm
SUT-457SP R

1.5 mm
SUT-458SP R

2.0 mm
SUT-459SP R

US 2-pin connector

Total length* 15.5 cm (6.1")

Working length* 5.7 cm (2.25")



0.4 mm
SUT-460SP R

0.7 mm
SUT-463SP R

1.0 mm
SUT-465SP R

1.2 mm
SUT-467SP R

1.5 mm
SUT-468SP R

2.0 mm
SUT-469SP R

US 2-pin connector

Total length* 19.1 cm (7.5")

Working length* 9.5 cm (3.75")



0.4 mm
SUT-530SP R

0.7 mm
SUT-533SP R

1.0 mm
SUT-535SP R

1.2 mm
SUT-537SP R

1.5 mm
SUT-538SP R

2.0 mm
SUT-539SP R

US 2-pin connector

Total length* 22.1 cm (8.7")

Working length* 12.7 cm (5.0")



0.4 mm
SUT-540SP R

0.7 mm
SUT-543SP R

1.0 mm
SUT-545SP R

1.2 mm
SUT-547SP R

1.5 mm
SUT-548SP R

2.0 mm
SUT-549SP R


US 2-pin connector

Total length* 24.6 cm (9.7")

Working length* 15.2 cm (6.0")

Sold through Sutter Exclusive Partners. See back cover for contact details.

* Lengths for orientation purposes; may vary slightly.



Disclaimer:

Product availability is subject to regulatory approval in individual markets.
Products may therefore not be available in all markets.



Manufactured for

SUTTER MEDICAL TECHNOLOGIES USA, INC.

120 INTERSTATE NORTH PARKWAY · SUITE 118 · ATLANTA, GA 30339

WWW.SUTTER-USA.COM · INFO@SUTTER-USA.COM