



Calvian[®]

Bipolar Forceps




PRECISION
ELECTROSURGERY

Calvian® Features

- ✓ Insulated shaft and handle
- ✓ Standard grasping forceps design
- ✓ Models with or without suction available
- ✓ Autoclavable
- ✓ 15° and 45° angled tips

15° and 45° angled tips



Calvian® duckbill+




Calvian® duckbill+



Calvian®




Calvian®



Calvian® duckbill⁺



Front-end coagulation is possible.



Hard-to-reach vessels and structures can be reached thanks to the recessed opening of the suction channel.

15° and 45° angled tips



Calvian® duckbill⁺



Calvian® duckbill⁺



The design of the tips enables front-end coagulation.



Hard-to-reach vessels and tissue can be reached without being grasped.



The position of the suction channel provides room for grasping of fine structures and helps to avoid unwanted suction.



Calvian® duckbill+

Calvian® duckbill+



12 cm (4.75")

Working length
12 cm (4.75")

15° angled tip
duckbill+
With suction

70 09 38



12 cm (4.75")

Working length
12 cm (4.75")

45° angled tip
duckbill+
With suction

70 09 39

Calvian® To Scale 1:1



Calvian® To Scale 1:1



18 cm (7.0")
Working length
18 cm (7.0")

15° angled tip
duckbill
Without suction

70 09 49



18 cm (7.0")
Working length
18 cm (7.0")

45° angled tip
duckbill
Without suction

70 09 48

Calvian® To Scale 1:1



23 cm (9.0")

Working length
23 cm (9.0")

45° angled tip
duckbill
With suction

70 09 46

23 cm (9.0")

Working length
23 cm (9.0")

45° angled tip
duckbill
Without suction

70 09 47

Accessories

Generator connector	Length	Bipolar cable 2 x 2 mm pins
Valleylab/Covidien	4,5 m	37 01 50 R
US Standard, ERBE ICC, International Bipolar cable 2 x 2 mm pins Length	4,5 m	37 01 35 R
Olympus, Söring, Coa Comp	4,5 m	37 01 53 R
Bipolar adaptor Sutter, Length: 40 cm For: Europ. flat connector		37 01 00 R



SUTTER MEDICAL TECHNOLOGIES USA, INC.

120 INTERSTATE NORTH PARKWAY · SUITE 118 · ATLANTA, GA 30339

WWW.SUTTER-USA.COM · E-MAIL: INFO@SUTTER-USA.COM

CE 0297