



Swyng[®] non-stick

Bipolar Forceps
single-use



non-stick Technology

Our proprietary non-stick technology was developed to reduce sticking and charring during coagulation – thanks to optimal heat transfer capabilities.

Low weight and optimal balance

Light-weight aluminum core and innovative product design featuring a dynamic curved shape that follows the natural ergonomics of the hand for improved comfort and ideal balance.



Two Product Lines



Standard

The pylon guide prevents scissoring of the tips and facilitates the precise coagulation and dissection in your surgical cases. Available with **straight** and **angled tips**.




Slim

The slim profile was designed for improved visualization of the surgical field, even in narrow operating corridors.

Standard

Large Product Portfolio


- 18 different patterns
- straight and angled tip styles
- three lengths and tip sizes available
- two product lines



Product availability is subject to regulatory approval in individual markets. Products may therefore not be available in all markets. The listed lengths and sizes serve as a guideline and may be rounded up or down. The actual lengths may vary slightly.




Slim




0.5 mm
78 44 105 SL

Slim



1.0 mm
78 44 110 SL

Slim




1.5 mm
78 44 115 SL

US 2-pin connector	
Total length	21.0 cm (8.25")
Working length	8.0 cm (3.0")




Slim




0.5 mm
78 44 405 SL

Slim



1.0 mm
78 44 410 SL

Slim



1.5 mm
78 44 415 SL

US 2-pin connector	
Total length	24.0 cm (9.5")
Working length	11.0 cm (4.25")



single-use



Standard

0.5 mm
78 44 905 ST

Standard

1.0 mm
78 44 910 ST

Standard

1.5 mm
78 44 915 ST

US 2-pin connector	
Total length	19.0 cm (7.5")
Working length	6.0 cm (2.5")



Standard

1.0 mm
45° angled up
78 44 710 ST

Standard

1.0 mm
45° angled down
78 44 810 ST

Standard

0.5 mm
78 44 105 ST

Standard

1.0 mm
78 44 110 ST

Standard

1.5 mm
78 44 115 ST

US 2-pin connector	
Total length	21.0 cm (8.25")
Working length	8.0 cm (3.0")



Standard

1.0 mm
45° angled up
78 44 610 ST

Standard

0.5 mm
78 44 405 ST

Standard

1.0 mm
78 44 410 ST

Standard

1.5 mm
78 44 415 ST

US 2-pin connector	
Total length	24.0 cm (9.5")
Working length	11.0 cm (4.25")




**DYNAMIC
CURVED
SHAPE**

Our disposable
non-stick bipolar forceps line!

The right instrument for your surgical needs

The Swyng® non-stick bipolar forceps offer precision and versatility in the operating room with a wide range of available models!



© Sutter Medizintechnik · Subject to change · REF 1329 B_USA – V10 · printed on acid free paper

SUTTER MEDICAL TECHNOLOGIES USA, INC.

120 INTERSTATE NORTH PARKWAY · SUITE 118 · ATLANTA, GA 30339

WWW.SUTTER-USA.COM · INFO@SUTTER-USA.COM