

# Calvian endo-pen<sup>®</sup>


**non-stick** bipolar forceps



Combining the benefits of the Calvian endo-pen<sup>®</sup> bipolar forceps with our established SuperGliss<sup>®</sup> **non-stick technology**

[Laboratory tests confirm the non-stick properties.\*]


**New!**



lateral view

### **non-stick technology**

Our specially developed material prevents overheating of the tips during coagulation.



Calvian endo-pen<sup>®</sup> bipolar forceps for tissue removal and bleeding control

- ✓ Slender shaft and handle
- ✓ No assembly, no loose parts
- ✓ Forceps-like handle
- ✓ Sturdy tips for coagulation and tissue removal

\* Data on file, Sutter Medizintechnik GmbH, Freiburg (Germany)

## Product Information

Instrument tray for one Calvian endo-pen® non-stick included.

### Calvian endo-pen® non-stick bipolar forceps

0.7 mm tip width  
15° angled tip  
working length: 10 cm (3.9")  
total length: 22 cm (8.6")

US pin connector  
70 09 89 SGS



horizontal movement



### Calvian endo-pen® non-stick bipolar forceps

0.7 mm tip width  
15° angled tip  
working length: 18 cm (7.0")  
total length: 30 cm (12.0")

US pin connector  
70 09 88 SGS



horizontal movement




## Accessories

Instrument trays available to hold one, two, five, or ten Calvian endo-pen® non-stick bipolar forceps.

### Instrument tray for one (1) Calvian endo-pen® non-stick bipolar forceps (without instruments)

70 17 78-01  
Connection for the rinsing port of the washer-disinfector included.




### Instrument tray for two (2) Calvian endo-pen® non-stick bipolar forceps (without instruments)

70 17 78-02  
Connection for the rinsing port of the washer-disinfector included.




### Instrument tray for five (5) Calvian endo-pen® non-stick bipolar forceps (without instruments)

70 17 78-05  
Connection for the rinsing port of the washer-disinfector included.



### Instrument tray for ten (10) Calvian endo-pen® non-stick bipolar forceps (without instruments)

70 17 78-10  
Connection for the rinsing port of the washer-disinfector included.



Cables for Calvian endo-pen® non-stick bipolar forceps can be found in the Calvian endo-pen® bipolar forceps product catalog.

Product availability is subject to regulatory approval in individual markets. Products may therefore not be available in all markets.



SUTTER MEDICAL TECHNOLOGIES USA, INC.

120 INTERSTATE NORTH PARKWAY · SUITE 118 · ATLANTA, GA 30339

WWW.SUTTER-USA.COM · INFO@SUTTER-USA.COM