




# Adson Bipolar Forceps



0.7 mm  
SUT-377NS

US 2-pin connector  
Total length\* 12.7 cm (5.0")

1.0 mm  
SUT-380NS

US 2-pin connector  
Total length\* 12.7 cm (5.0")

1.5 mm  
SUT-382NS

US 2-pin connector  
Total length\* 12.7 cm (5.0")



## non-stick technology

Silver alloy tips prevent overheating of the tips during coagulation. The non-stick effect is permanently ensured, as the tips are made of solid noble metals.




## Handle

The ergonomically designed handle provides less fatigue and a well balanced instrument.

Sold through Sutter Exclusive Partners. See back cover for contact details.

\* Lengths for orientation purposes; may vary slightly.



Disclaimer:

Product availability is subject to regulatory approval in individual markets.  
Products may therefore not be available in all markets.

**Manufactured for**

**SUTTER MEDICAL TECHNOLOGIES USA, INC.**

120 INTERSTATE NORTH PARKWAY · SUITE 118 · ATLANTA, GA 30339

WWW.SUTTER-USA.COM · INFO@SUTTER-USA.COM