



Bipolar Irrigation Forceps



Bipolar irrigation forceps · K-YASARGIL non-stick · hidden irrigation



hidden irrigation



non-stick technology

Silver alloy tips prevent overheating of the tips during coagulation. The non-stick effect is permanently ensured, as the tips are made of solid noble metals.



Hidden irrigation

Developed to increase line of sight compared to regular irrigation. The hidden irrigation allows maximum visibility.



Irrigation cap

A universal hub connector allows irrigation tubing to be connected easily.



Handle

The ergonomically designed handle provides less fatigue and a well balanced instrument.



0.5 mm
SUT-481NS HR



0.7 mm
SUT-483NS HR



1.0 mm
SUT-485NS HR



1.3 mm
SUT-488NS HR

US 2-pin connector

Total length* 15.9 cm (6.25")

Working length* 6.0 cm (2.36")



0.5 mm
SUT-491NS HR



0.7 mm
SUT-493NS HR



1.0 mm
SUT-495NS HR



1.3 mm
SUT-498NS HR

US 2-pin connector

Total length* 18.1 cm (7.12")

Working length* 8.3 cm (3.25")



- 0.5 mm
SUT-501NS HR
- 0.7 mm
SUT-503NS HR
- 1.0 mm
SUT-505NS HR
- 1.3 mm
SUT-508NS HR

US 2-pin connector

Total length* 20.0 cm (7.88")

Working length* 10.2 cm (4.0")




- 0.5 mm
SUT-521NS HR
- 0.7 mm
SUT-523NS HR
- 1.0 mm
SUT-525NS HR
- 1.3 mm
SUT-528NS HR

US 2-pin connector

Total length* 27.9 cm (11.0")

Working length* 17.2 cm (6.75")

* Lengths for orientation purposes; may vary slightly.



Disclaimer:

Product availability is subject to regulatory approval in individual markets.
Products may therefore not be available in all markets.

© Sutter Medizintechnik · Subject to change · REF 1278_USA – T02 · printed on acid free paper



Manufactured for
SUTTER MEDICAL TECHNOLOGIES USA, INC.
120 INTERSTATE NORTH PARKWAY · SUITE 118 · ATLANTA, GA 30339
WWW.SUTTER-USA.COM · INFO@SUTTER-USA.COM