




Calvian endo-pen[®]

Bipolar Forceps



PRECISION
ELECTROSURGERY



Horizontal movement

To Scale 1:1



US pin connector	
Total length*	22.0 cm (8.6")
Working length*	10.0 cm (3.9")

* Lengths for orientation purposes; may vary slightly.

15°
0.7 mm tip width
15° angled tip
US pin connector
70 09 59 S
lateral view


15°
1.2 mm tip width
15° angled tip
US pin connector
70 09 67 S
lateral view

15°
1.5 mm tip width
15° angled tip
US pin connector
70 09 69 S
lateral view

45°
0.7 mm tip width
45° angled tip
US pin connector
70 09 92 S
lateral view

45°
1.2 mm tip width
45° angled tip
US pin connector
70 09 94 S
lateral view

45°
2.5 mm tip width
45° angled tip
US pin connector
70 09 55 S
lateral view



US pin connector	
Total length*	30.0 cm (12.0")
Working length*	18.0 cm (7.0")

* Lengths for orientation purposes; may vary slightly.

15°
0.7 mm tip width
15° angled tip
US pin connector
70 09 86 S
lateral view


15°
1.2 mm tip width
15° angled tip
US pin connector
70 09 87 S
lateral view

15°
1.5 mm tip width
15° angled tip
US pin connector
70 09 88 S
lateral view

45°
0.7 mm tip width
45° angled tip
US pin connector
70 09 89 S
lateral view

45°
1.2 mm tip width
45° angled tip
US pin connector
70 09 90 S
lateral view

45°
2.5 mm tip width
45° angled tip
US pin connector
70 09 56 S
lateral view



Vertical movement

To Scale 1:1



US pin connector	
Total length*	22.0 cm (8.6")
Working length*	10.0 cm (3.9")

* Lengths for orientation purposes; may vary slightly.



top view

15°
1.2 mm tip width
15° angled tip
US pin connector
70 09 68 S




top view

45°
2.5 mm tip width
45° angled tip
US pin connector
70 09 83 S



top view

70°
2.5 mm tip width
70° angled tip
US pin connector
70 09 84 S



US pin connector	
Total length*	30.0 cm (12.0")
Working length*	18.0 cm (7.0")

* Lengths for orientation purposes; may vary slightly.



top view

15°
0.7 mm tip width
15° angled tip
US pin connector
70 09 57 S



top view

15°
1.2 mm tip width
15° angled tip
US pin connector
70 09 58 S



top view

15°
2.5 mm tip width
15° angled tip
US pin connector
70 09 52 S



top view

45°
2.5 mm tip width
45° angled tip
US pin connector
70 09 50 S




top view


70°
2.5 mm tip width
70° angled tip
US pin connector
70 09 54 S

Advantages


- ✓ Forceps-like handle
- ✓ No assembly, no loose parts
- ✓ Slender shaft and handle
- ✓ Tip widths from 0.7 mm to 2.5 mm available
- ✓ Different angles: 15°, 45°, and 70°
- ✓ Reusable and autoclavable
- ✓ Each Calvian endo-pen® comes with an instrument tray for safe cleaning and sterilization.



Calvian endo-pen® bipolar forceps for tissue removal and bleeding control.








Forceps jaws




Paddle-style jaws

Accessories

Generator connector	Length	Safety connector / US pin connector
	3.0 m	37 01 54 S
	4.5 m	37 01 35 S
	4.5 m	37 01 50 S
	4.5 m	37 01 53 S
	3.0 m	37 01 40 S

Instrument tray, holds 1 Calvian endo-pen® (without instruments)

Connection for the rinsing port of the washer-disinfector included



70 17 78-01 for 1 instrument

Instrument tray, holds 2 Calvian endo-pen® (without instruments)


Connection for the rinsing port of the washer-disinfector included



70 17 78-02 for 2 forceps

Instrument tray, holds 5 Calvian endo-pen® (without instruments)


Connection for the rinsing port of the washer-disinfector included



70 17 78-05 for 5 forceps


Instrument tray, holds 10 Calvian endo-pen® (without instruments)

Connection for the rinsing port of the washer-disinfector included



70 17 78-10 for 10 forceps

Disclaimer: Products shown in this catalog are available in the US. Availability in other markets may vary.



© Sutter Medizintechnik · Subject to change · REF 1228 D_USA – R11 · printed on acid free paper



SUTTER MEDICAL TECHNOLOGIES USA, INC.
120 INTERSTATE NORTH PARKWAY · SUITE 118 · ATLANTA, GA 30339
WWW.SUTTER-USA.COM · INFO@SUTTER-USA.COM